

TECHSPEC[®] Aperture Kit for Achromat Prototyping Kit



Aperture Kit for Achromat Prototyping Kit

Stock #54-720 3-5 DAYS



1



A\$233.76

ADD TO CART

Qty 1-5

A\$233.76

Qty 6+

A\$217.44

Volume Pricing

[Request Quote](#)

Get Product Downloads

Type:
Prototyping Kit

Function:
Kit

Regulatory Compliance

RoHS:

Compliant

Certificate of Conformance:
[View](#)
