

TECHSPEC® C-Mount Male to Nikon (M25) Female Step-Down Adapter



C-Mount to Nikon (M25) Adapter (#11-148)

Stock **#11-148** **20+ In Stock**

A\$83²⁰

ADD TO CART

Volume Pricing	
Qty 1-4	A\$83.20 each
Qty 5-9	A\$74.88 each
Qty 10-25	A\$70.72 each
Need More?	Request Quote

Product Downloads

General

Objective Mount **Type:**

Objective **Type of Optics:**

Objective Adapter **Function:**

Physical & Mechanical Properties

