

TECHSPEC®

Cell with Two Retainers for Achromat Prototyping Kit



Cell with Two Retainers for Achromat Prototyping Kit

Stock #58-255 3-5 DAYS



1



A\$233.76

ADD TO CART

Qty 1-5

A\$233.76

Qty 6+

A\$205.92

Volume Pricing

[Request Quote](#)

Get Product Downloads

Type:
Prototyping Kit

Function:
Kit

Regulatory Compliance

RoHS:

Compliant

Certificate of Conformance:
[View](#)
