

TECHSPEC[®] FITC Filter Cube Set - Zeiss



Olympus Filter Cube #86-833



Stock **#67-004-ZSS** NEW [CONTACT US](#)

-

1

+

A\$2,160⁰⁰

ADD TO CART

Volume Pricing	
Qty 1-9	A\$2,160.00 each
Qty 10+	A\$1,944.00 each
Need More?	Request Quote

Product Downloads

General	
Fluorescence Filter Kit	Type:
FITC	Compatible Fluorophore:
#67-028	Excitation Filter:

#67-031	Emission Filter:
#67-080	Dichroic Filter:
Zeiss and Zaber	Compatible Microscope Brand:
Optical Properties	
Fused Silica (Corning 7980)	Substrate: □
Hard Coated	Coating:
467 - 498	Excitation Wavelength (nm):
513 - 556	Emission Wavelength (nm):
506.00	Dichroic Cut-On Wavelength (nm):
467 - 556	Wavelength Range (nm):
Regulatory Compliance	
View	Certificate of Conformance:

Product Details

- Ensure Proper Alignment of Filter Sets
- Easily Adapt to Common Miroscopes
- Available in a Variety of Wavelengths
- Unmounted [TECHSPEC® Fluorescence Filter Kits](#) Also Available

TECHSPEC® Pre-Mounted Fluorescence Filter Cube Sets are designed to easily integrate into popular Olympus, Nikon, Zeiss, and Zaber microscopes. Each fluorescence filter cube set consists of one [TECHSPEC Fluorescence Dichroic Filter](#) and two [TECHSPEC Fluorescence Bandpass Filters](#) that have been specifically oriented to achieve optimum excitation and emission of popular fluorophores. Additionally these filters are also used in polymerase chain reaction (PCR) diagnostic instruments. For increased adaptability, filters can be easily removed and replaced to meet the needs of the application. Unmounted [TECHSPEC Fluorescence Filter Kits](#) are also available.

Built-to-Order Product: Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).

Built-to-Order Product: Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).