

TECHSPEC® Qdot® 625 Filter Cube Set - Zeiss



Olympus Filter Cube #86-833



Stock #86-376-ZSS NEW CONTACT US

-

1

+

A\$2,160⁰⁰

ADD TO CART

Volume Pricing	
Qty 1-9	A\$2,160.00 each
Qty 10+	A\$1,944.00 each
Need More?	Request Quote

Product Downloads

General	
Fluorescence Filter Kit	Type:
Qdot® 625	Compatible Fluorophore:
#86-349	Excitation Filter:

#86-357	Emission Filter:
#86-332	Dichroic Filter:
Zeiss and Zaber	Compatible Microscope Brand:
Optical Properties	
Fused Silica (Corning 7980)	Substrate: □
Hard Coated	Coating:
415 - 455	Excitation Wavelength (nm):
618 - 633	Emission Wavelength (nm):
510.00	Dichroic Cut-On Wavelength (nm):
415 - 633	Wavelength Range (nm):
Regulatory Compliance	
View	Certificate of Conformance:

Product Details

- Ensure Proper Alignment of Filter Sets
- Easily Adapt to Common Microscopes
- Available in a Variety of Wavelengths
- Unmounted [TECHSPEC® Fluorescence Filter Kits](#) Also Available

TECHSPEC® Pre-Mounted Fluorescence Filter Cube Sets are designed to easily integrate into popular Olympus, Nikon, Zeiss, and Zaber microscopes. Each fluorescence filter cube set consists of one [TECHSPEC Fluorescence Dichroic Filter](#) and two [TECHSPEC Fluorescence Bandpass Filters](#) that have been specifically oriented to achieve optimum excitation and emission of popular fluorophores. Additionally these filters are also used in polymerase chain reaction (PCR) diagnostic instruments. For increased adaptability, filters can be easily removed and replaced to meet the needs of the application. Unmounted [TECHSPEC Fluorescence Filter Kits](#) are also available.

Built-to-Order Product: Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).

Built-to-Order Product: Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).