

## TECHSPEC® TRITC Filter Cube Set - Olympus



Olympus Filter Cube #86-833 with Pre-Mounted Fluorescence Filters



Stock #67-007-OLY 3-5 DAYS

⊖ 1 ⊕ **A\$1,856.00**

**ADD TO CART**

Qty 1-9

**A\$1,856.00**

Qty 10+

**A\$1,688.00**

Volume Pricing

[Request Quote](#)

Get Product Downloads

**Substrate:**  [Fused Silica](#) (Coming 7980)

**Coating:**  
Hard Coated

**Type:**  
Fluorescence Filter Kit

**Excitation Wavelength (nm):**  
532 - 554

**Compatible Fluorophore:**  
TRITC

**Emission Wavelength (nm):**  
570 - 613

**Dichroic Cut-On Wavelength (nm):**  
562.00

**Excitation Filter:**  
[#67-032](#)

**Emission Filter:**  
[#67-034](#)

**Dichroic Filter:**  
[#67-082](#)

**Wavelength Range (nm):**  
532 - 613

**Compatible Microscope Brand:**  
Olympus

## Regulatory Compliance

**RoHS:**  
Not Compliant

**Certificate of Conformance:**  
[View](#)

---